

## Upgrade POLYVAR MET / - '2/ - `SC / MeF3

### Video/Camera-Adapter

C-Mount/A30 vertical, inclusive optic for CCD ½"

- Port vertically instead of Polaroid case, rotate able 360°

C-Mount/A40 vertical, inclusive optic for CCD • "

- Port vertically instead of Polaroid case, rotate able 360°

B-Mount/A50 vertical, inclusive optic for Sony bayonet

- Port vertically instead of Polaroid case, rotate able 360°

C-Mount/A50 lateral inclusive optic also for Reichert MeF3 microscopes

- Port vertically instead of 35mm camera, rotate able 360°  
(At MeF3 probably have to screw off/rework cover on front/left!)

C-Mount/A120 vertical, inclusive optic and cover

- for CCD camera on Polyvar Met without Polaroid port  
- Port vertically instead of standard cover, rotate able 360°

### Real View Module for Incident Light

The Real View Optic corrects the lateral and vertically fault. Engaged there is a shift of 30 mm to the objective corpus and additional we receive a 6 inch out of a 4 inch scope and/or 8 inch a 6 inch microscope.

At the magnification wheel 0,8x the field of view reduces from 24 to 20mm.

Stage modification and transmitted light on require.

### Motorized Revolving Nosepiece

Modification of manual nosepiece to a motorized nosepiece for particle free and ergonomic operation, including two direction buttons, circuit, power supply, motor and gearing with holder, 6 objective stops and special installation.

